



ACRO CLUBE DA MAIA

X Maia International Acro Cup

4 - 6 mar, 2016, Maia (POR)

Qualifying
General Results



Page. 1

Mixed Pair/Senior

Rk		Time	EXE	ART	DIFF	PEN	Score	Total	GAP
1	SEA  Bruno Tavares/ Liana Asseiceiro								
	Balance	(2:30)	7,450	8,350	10,030 (153)		25,830	52,150	
	Dynamic	(2:26)	8,200	8,200	9,920 (136)		26,320		
2	APG  Marieke Uytterhoeven/ Maxime Vantournhout								
	Balance	(2:28)	8,100	7,950	9,600 (120)		25,650	51,590	0,560
	Dynamic	(2:28)	8,100	8,400	9,440 (112)		25,940		
3	AAC  M ^a Côdeas/ Pedro Calheiros								
	Balance	(2:31)	7,850	7,900	10,000 (150)	0,3	25,450	50,020	2,130
	Dynamic	(2:26)	7,000	7,750	9,820 (131)		24,570		
4	TSGTSV  Melissa Maier/ Patrick Schönholz								
	Balance	(2:26)	7,200	7,300	9,960 (146)	0,9	23,560	48,330	3,820
	Dynamic	(2:06)	7,700	7,450	9,620 (121)		24,770		
5	ARTG  Jasper Criel/ Lynn Eronini								
	Balance	(2:30)	7,400	7,600	9,600 (121)	0,3	24,300	47,150	5,000
	Dynamic	(2:24)	7,800	7,550	7,500 (71)		22,850		
6	SKYL  Levi Gooris/ Mirte Looock								
	Balance	(2:31)	7,100	7,450	9,540 (119)	0,6	23,490	47,120	5,030
	Dynamic	(2:22)	7,000	7,350	9,280 (107)		23,630		
7	APOMB  Bruno Felizardo/ Camila Silva								
	Balance	(2:21)	7,700	7,000	8,500 (91)	0,3	22,900	45,300	6,850
	Dynamic	(2:04)	7,350	7,150	7,900 (78)		22,400		
8	DINAM  Guillem Martinez/ Silvia Coronado								
	Balance	(2:23)	6,550	7,100	8,200 (85)	2,2	19,650	45,200	6,950
	Dynamic	(2:16)	8,450	8,400	8,700 (95)		25,550		
9	ACRBT  Adrian Hernandez/ Carolina Meija								
	Balance	(2:30)	6,800	7,250	8,400 (89)	1,2	21,250	43,400	8,750
	Dynamic	(2:31)	6,850	7,200	8,100 (82)		22,150		



ACRO CLUBE DA MAIA

X Maia International Acro Cup

4 - 6 mar, 2016, Maia (POR)

Qualifying
General Results



Page. 2

Mixed Pair/Senior

10	NDGA		Emily Arnold/ George Lee						
			Balance	(2:29)	8,200	8,000	8,500 (91)	0,3	24,400 24,400 27,750
			Dynamic	(0:00)	0,000	0,000	0,000		0,000
	APG		Robbe Steenssens/ Yoni Berteloot						
			Balance	(0:00)	0,000	0,000	0,000		0,000 0,000 52,150
			Dynamic	(0:00)	0,000	0,000	0,000		0,000